

# TURBISCAN OIL SERIES



ASTM  
D7061



#### HEAVY FUEL STABILITY

ASTM D7061 compliant to determine the stability reserve of crude and heavy fuel oil in one click



#### MULTI APPLICATIONS

Stability, asphaltene dispersant efficiency, water in oil demulsification, compatibility of blends



#### PORTABLE AND ROBUST

Designed for R&D labs and on-field analysis (<5kg)



#### THE REFERENCE

Turbiscan technology is the worldwide reference for the colloidal stability analysis

## STABILITY ANALYZER FOR CRUDE & FUEL OILS

STABILITY & SIZE

[www.formulation.com](http://www.formulation.com)




## TURBISCAN®, THE REFERENCE

Turbiscan® is used world-wide to detect at an early stage all kinds of destabilization such as coalescence, flocculation, creaming, sedimentation, etc... Turbiscan® Oil Series allows the stability analysis of heavy fuel oil thanks to an innovative method granted ASTM D7061.



### MEASUREMENT PRINCIPLE

2000+ publications  
700+ patents  
50+ countries




Turbiscan® works on Static Multiple Light Scattering in both Transmission (T) and Backscattering (BS), in order to analyze low and high concentration dispersions. T & BS signals depend on particle size and concentration: **BS & T = f ( d / φ )**

The Turbiscan® OIL SERIES acquires T & BS every 20 microns along the sample height. Scans are repeated during ageing time to detect any variation of the signal due to a destabilization, such as particle migration and/or particle size variation.

### D7061 ASTM METHOD STABILITY OF HEAVY FUEL OILS

Quick and easy method to determine the stability reserve upon ageing simulation (precipitation of asphaltene in alkane)



### KEY BENEFITS

#### ACCURATE STABILITY ANALYSIS

- Measurement of stability reserve thanks to ageing simulation (precipitation of asphaltene in alkane)
- High repeatability, accuracy and objectivity

#### FAST STABILITY ANALYSIS

- Up to 200 times faster than naked eye
- Only 15 minutes to determine stability of heavy fuels and crude oils
- One product, One stability, One number

#### EASY STABILITY ANALYSIS

- State of the art software with 1-click data treatment and report generation
- No calibration, no parameters, simply dilute in toluene and precipitate with heptane
- Robust and portable device to carry on the field if necessary

### APPLICATIONS



Oil & Petroleum



Heavy fuel



Crude Oils



Effects of additives



Effects of blending

### TECHNICAL SPECIFICATIONS

Cell Volume	7 ml
ASTM D-7061 compliant	•
ISO TR 13097 compliant	•
Quantitative monitoring of dispersion stability	•
Migration velocity & hydrodynamic diameter	•
Size range	10 nm to 1 mm
Concentration range	0.0001 to 60% v/v
Repeatability auto	± 0.1 %
Repeatability manual	± 0.25 %
Temperature control	No
Dimensions (cm)	34x16x29
Weight (kg)	5

